

Publications

Study

Users report

James R. Dunn, Carlos A. Munoz  
Loma Linda University, School of Dentistry,  
Center of Dental Research, Loma Linda/USA

## Clinical evaluation of a new microhybrid composite (Venus®)

# Venus®

## Scientific Information

# James R. Dunn Carlos A. Munoz

Loma Linda University, USA

## Objective

The purpose of this clinical trial was to evaluate the *in vivo* performance of a new generation of composite resin (micro-hybrid technology), Venus<sup>®</sup>, for anterior restorations over a one-year period.


## Materials and Methods

A total of 53 restorations (classes III, IV, V, diastema restorations, veneers) using Venus<sup>®</sup> were placed in the anterior teeth of 28 subjects by the Center of Dental Research at the Loma Linda University School of Dentistry. Two dentists placed the restorations over a 60-day period. Gluma<sup>®</sup> Comfort Bond + Desensitizer was used as the dental adhesive.

At one-year, 46 restorations were evaluated. They were evaluated by one evaluator (dentist) for anatomic form, color match, marginal adaptation, retention, marginal discoloration, surface staining, secondary caries and polishability by using modified USPHS evaluation criteria.

## Results

	Alpha	Bravo	Charlie	Delta	Eta
Anatomic Form (wear)	100%	0%	0%	0%	0%
Color Match	93%	7%	0%	0%	0%
Marginal Adaptation	100%	0%	0%	0%	0%
Retention	100%	0%	0%	0%	0%
Marginal Discoloration	100%	0%	0%	0%	0%
Surface Staining	100%	0%	0%	0%	0%
Secondary Caries	100%	0%	0%	0%	0%
Polishability	54%	46%	0%	0%	0%



## Conclusions

In the categories anatomic form, marginal adaptation, retention, marginal discoloration, surface staining and secondary caries all restorations were ranked Alpha.

None of the subjects exhibited post-operative sensitivity or exaggerated gingival response.

Based on the one-year findings of this study, the Venus<sup>®</sup> micro-hybrid composite is an excellent choice for anterior restorations.

### Contact in Germany

**Heraeus Kulzer GmbH**  
Grüner Weg 11  
63450 Hanau  
Phone: +49 (0) 6181 355 444  
Fax: +49 (0) 6181 353 461  
info.dent@heraeus.com  
www.heraeus-kulzer.de

### Contact in the United Kingdom

**Heraeus Kulzer Ltd.**  
Heraeus House, Albert Road  
Northbrook Street, Newbury  
Berkshire, RG14 1DL  
Phone: +44 1 635 30500  
Fax: +44 1 635 30606  
Mail: sales@kulzer.com.uk  
www.heraeus-kulzer.com

### Contact in Australia

**Heraeus Kulzer Australia Pty. Ltd.**  
Locked Bag 750  
Roseville NSW 2069  
Phone: +61 29 417 8411  
Fax: +61 29 417 5093  
Mail: sales@kulzer.com.au  
www.kulzer.com.au