# Venus<sup>®</sup> *flow*

Ⓚ	Gebrauchsanweisung .....	3
Ⓜ	Instructions for use .....	9
Ⓕ	Mode d'emploi .....	15
Ⓔ	Instrucciones de uso .....	21
Ⓘ	Istruzioni per l'uso .....	27
Ⓟ	Instruções de uso .....	33
Ⓝ	Gebruiksaanwijzing .....	39
Ⓢ	Bruksanvisning .....	45
Ⓓ	Brugsvejledning .....	51
Ⓝ	Bruksanvisning .....	57
Ⓕ	Käyttöohje .....	63
Ⓖ	Οδηγίες Χρήσης .....	69

**Heraeus**

# Venus<sup>®</sup> flow

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# Venus<sup>®</sup> flow

Instructions  
for use

Venus Flow is a flowable, light curing, radiopaque hybrid composite used for adhesive, tooth-coloured anterior and posterior restorations. Due to its low viscosity Venus Flow provides excellent coverage of the tooth structure.

## **Composition:**

Venus Flow is based on multifunctional (meth-) acrylate monomers and contains approx. 62% inorganic fillers (Ba-Al-F glass, SiO<sub>2</sub>) by weight.

## **Shades:**

A1, A2, A3, A3.5, A4, B2, B3, OA2, SB1, SB2, SBO, T2, HKA2.5, Baseline white

## **Applications:**

- Fissure sealing
- Enlarged fissure sealing
- Class V fillings
- Minimally invasive Class I and II fillings in areas not subjected to masticatory forces
- Minimally invasive Class III fillings
- Smaller contour and shade adjustment on the enamel and dentine

- Small surface repairs to direct and indirect restorations in combination with a suitable bonding agent
- Splinting of teeth
- Cavity lining
- Bracket retention

### **Kind of application:**

Prepare the tooth structure for treatment by cleaning it with a fluoride free polishing paste. Before treatment, while the teeth are still moist, select the Venus Flow shade using the Venus 2Layer shade guide.

1. Remove any caries and complete minimally invasive preparation of the tooth structure according to the general guidelines for adhesive fillings. Then clean and dry the cavity. The use of a rubber dam is recommended.
2. A quality adhesive system (e.g. iBond®, GLUMA® Comfort® Bond + Desensitizer, GLUMA® Solid Bond) should be applied to the enamel and dentin in accordance with the manufacturer's instructions before inserting the filling material. Venus Flow can be applied without adhesive as fissure sealing to purely enamel areas directly following enamel etching (e.g. with GLUMA® Etch 20, GLUMA® Etch 35). Avoid any contamination of conditioned or adhesive-treated surfaces with moisture, blood or saliva.

3. Apply Venus Flow in thin layers (max. 2 mm, baseliner max. 1 mm) to the cavity. When applying Venus Flow ensure that it flows slowly and evenly to prevent air bubbles – apply a gentle stream of air to facilitate maximum marginal contact. If necessary, remove any air bubbles with a probe before polymerisation.
4. Polymerise each layer according to the table below using a suitable light-curing unit (e.g. Translux<sup>®</sup>) for 20 or 40 sec. A dispersion layer forms during polymerisation. This layer should not be touched or removed if other composite layers are to be applied.

Shades	Time
A1, A2, A3, A3.5, B2, B3, SB1, SB2, T2, HKA2.5	20 s
A4, OA2, SBO, Baseliner	40 s

5. Venus Flow can be prepared and polished immediately following polymerisation. Finishing diamonds, tungsten carbide instruments, flexible discs, silicone polishers and polishing brushes are suitable for preparing.

During long restorative procedures, the surgery light should be temporarily directed away from the treatment area to prevent premature initiation of the polymerisation reaction.

**Note:**

- A Heraeus Kulzer Translux® light unit or a light unit with comparable intensity (min. 400-500 mW/cm<sup>2</sup>) should be used. The same polymerisation time is required with an LED light unit. If a stronger polymerisation unit (e.g. a plasma light unit with a greater output than 1200 mW/cm<sup>2</sup>) is used, the polymerisation time can be reduced by 25%.
- The polymerisation times given in the instructions for use are the minimum times required to attain optimum polymerisation.
- During polymerisation the light guide should be placed as close as possible to the surface being polymerised.
- The tips supplied for precise, hygienic application are intended for single-use and should be discarded after use. After removing the material from the syringe take the pressure off the plunger to prevent continued flow. Replace the cap of the syringe after use.
- With small surfaces in caries-free exposed pulp cavities, the exposed area should be directly capped before beginning the adhesive treatment (e.g. a combination of CaOH and glass ionomer cement). As small an area as possible of the dentine near the pulp should be covered with a cavity lining (e.g. glass ionomer cement).
- Do not use Venus Flow with eugenol cements, as these could inhibit polymerisation.

Avoid contact with the skin or eyes. On contact with the eyes, rinse with plenty of water and consult a doctor if necessary.

Do not use Venus Flow after the expiry date. Do not store above 25°C/77°F. Avoid direct sunlight.

Venus Flow should not be used with patients who have an allergic reaction to any of its constituents.

### **Delivery units / Packing units**

Syringes 1 x 1,8 g

Assortments as syringes 4 x 1,8 g

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to their suitability for the intended processes and uses. The application, use and processing of the products are beyond our control, and therefore, entirely your own responsibility. Should, in spite of this, liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery.

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Heraeus Kulzer GmbH
Heraeus Kulzer, Inc.