

Heraeus



Venus[®] Supra

Long-lasting, High Gloss Polishers

Venus[®]

The new aesthetics.



Feel the difference.
See the difference.

Venus Supra makes the difference.



Polished composite restoration



Venus[®] Supra
polished composite restoration



Venus Supra – The finishing touch for your restorative masterpiece.

New universal, two-step, diamond polishing system achieves long-lasting, high gloss polishing for every composite restoration – in only two steps for a naturally beautiful shine.

Venus Supra – Practicable and economical

Venus Supra polishers are autoclavable and can be reused several times depending on form and extension of the polished restorations.



Venus Supra – Polishers for “Supra” results

Use Venus Supra for polishing all types of composites – you will achieve long-lasting, high gloss polishing effects for naturally beautiful composite restorations!

Venus Supra pre- and high gloss polishers are made of silicone highly filled with microfine diamond powder. The five different shapes provide a universal solution for all anterior and posterior restorations.

Venus Supra – Two steps to long-lasting aesthetics

1. Red → Pre-Polishers



2. Grey → High gloss Polishers




Create invisible, naturally glossy composite restorations in only two steps without several changes of instruments and with low substance removal.

Venus Supra – Benefits at a glance

Highly filled with microfine diamond powder	Long-Lasting, High Gloss
Universally applicable	One system for all composite and restoration types
Only two steps	Time saving
Autoclavable and re-usable	Economical

Naturally beautiful restorations every time!

Venus Supra – High gloss for every composite



Source: Heraeus Kulzer R&D, Wehrheim, Germany
Unpublished. Data on file.

Venus Supra achieved the lowest surface roughness reducing plaque accumulation and pigmentation over time. The high gloss (reflection) showed by every composite polished with Venus Supra leads to long-lasting aesthetic results.



Do you want to know more about Venus?
Visit www.heraeus-venus.com

Venus Supra

Article _____ Order No. _____


Venus Supra Kit _____ 66039110
 10 polishers: each of all five shapes as pre- and high gloss polishers



Venus Supra Refills

10 polishers: one shape as pre- and high gloss polishers (each 5x)

■ Cup _____ 66039161



■ Disk _____ 66039162



■ Flame _____ 66039163



■ Large Cup _____ 66039164



■ Large Flame _____ 66039165



Contact in Germany
Heraeus Kulzer GmbH
 Grüner Weg 11
 63450 Hanau
 info.lab@heraeus.com
 www.heraeus-kulzer.com

Contact in the United Kingdom
Heraeus Kulzer Ltd.
 Heraeus House, Albert Road
 Northbrook Street, Newbury
 Berkshire, RG14 1DL
 Phone ++44 (0) 1635 30500
 Fax ++44 (0) 1635 30606
 sales.uk@heraeus.com
 www.heraeus-kulzer.com

Contact in Australia
Heraeus Kulzer Australia Pty. Ltd.
 Locked Bag 750
 Roseville NSW 2069
 Phone +61 29.417-8411
 Fax +61 29.417-5093
 sales@heraeus.com.au
 www.heraeus-kulzer.com